**CENTER OF EXPERTISE** 

### **ADVANCED X-RAY TOPOGRAPHY**

### INDUSTRY READY DETECTION OF TSDs AND BPDs IN SIC WAFERS

C. Reimann<sup>1</sup>, C. Kranert<sup>2</sup>, S. Kobayashi<sup>3</sup>, Y. Ueji<sup>3</sup>, K. Shimamoto<sup>3</sup>, K. Omote<sup>3</sup>, A. Vigliante<sup>1</sup>

<sup>1</sup>Rigaku Europe SE, Neu-Isenburg, Germany
<sup>2</sup>Fraunhofer IISB, Erlangen, Germany
<sup>3</sup>Rigaku Corp, Tokyo, Japan





Rigaku

000000

# accelerated representation map animated





Center of Expertise for X-Ray-Topography

### founded in 2021

- joint undertaking of **Rigaku Corporation** and **Fraunhofer IISB** to develop:
  - X-ray topography tools
  - measurement procedures
- defect recognition and counting algorithms
   for production, quality assurance and R&D application
- supporting from first demo measurements, tool training etc.



### Taking away the power from silicon...

- ... means leaving behind dislocation-free material!
- Dislocations can have an impact on device yield and reliability
- Every substrate is different



#### Industry requires semiconductor ready, reliable, high-throughput, non-destructive substrate characterization



### **Defects in silicon carbide**

Threading Screw Dislocations (TSD)



- Increase leakage current of diodes
- Detrimental to gate oxide lifetime
- Can result in growth pits causing problems during device processing

### Basal Plane Dislocations (BPD)



- Result in stacking faults under stress
- Reduce lifetime of
  - MOSFETs
  - Bipolar devices



# 211111 000-9455WEO Lansanioni **TSD Detection** X 0

### **TSD detection – "Old" vs. "New" language**

- TSD detection using chemical etching highly challenging
  - TED and TSD cannot be distinguished for standard KOH etching
  - Other alternative approaches come with different problems
- Use of (0008) reflex with XRTmicron
  - **TSD/TMD** appear as small, dark spots
  - Almost no other features visible (on today's production grade SiC)
  - Simple to classify and count → enables reliable and robust dislocation counting
    - Measurement reproducibility: 3%
    - Tool2Tool: 5% (depending on samples could be 10%)





### **TSD detection** Reproducibility

- Measurement of three neighboring wafers on both sides (Si-face and C-face)
- Differences between individual wafers can be resolved
- Differences between front and back side of wafer can be resolved
- From trend between data points, the distance between wafers can be estimated!







### **TSD detection** Reproducibility

- Full wafer measurement of many wafers from two SiC crystals
- Plot of TSD density over wafer number
- Monotonic TSD density variation is easily reproduced
- Even small differences (1%) between wafers can be resolved







### **TSD detection** Measurement times

- Industrial application:
  - Perfect accuracy usually not required
  - High throughput essential
- TSD detection requires minimum image quality → Throughput on full wafer scale is limited (A)
- Use of local measurement to increase throughput
  - Grid measurement (B) → only very small error at 3x throughput
  - Stripe measurement (C) → critically increased error at 4x throughput
    - Use of radial weighted average reduces error and narrowness of the distribution





### **TSD detection** SEMI Standard M91



### SEMI M91 - Test Method for Determination of Threading Screw Dislocation Density in 4H-SIC by X-Ray Topography

Volume(s): Materials Language: English Type: Single Standards Download (.pdf) SEMI Standards Copyright Policy/License Agreements

Member Price: \$113.00

Non-Member Price: \$150.00

Revision

### **READY FOR INDUSTRY!**



### **BPD Quantification**

0

DOC VESIVE

100000

.....

XXXX

### **BPD** quantification

Principle

- BPDs give only weak contrast in reflection geometry due to in-plane Burgers vector
  - → usually measured by XRT in transmission geometry
- Which properties to investigate?
  - BPD count?
  - BPD length?
- Volumetric density (cm/cm<sup>3</sup>) differs strongly (factor of 14 at 4° offcut) from etch pit density (cm<sup>-2</sup>)





### **BPD** quantification – "Noodle salad challenge" Principle

- Linear appearance of BPD pose challenges to counting
  - Detection of start and end of dislocation already challenging for low densities
  - At higher densities  $\rightarrow$  strong overlapping hinders differentiation between dislocations
  - Discrete counting (as for TSD) not possible



high density

low



### **BPD quantification – "Noodle salad challenge"** Principle

- Use of a calibrated approach to obtain results consistent to KOH results
  - Divide topogram into tiles (typically 1x1mm<sup>2</sup>)
  - Analyze topogram to obtain a numeric value
  - Etch wafer by KOH and obtain etch pit density for the same tiles
  - Calibration curve obtained by correlating both data sets





### **BPD** quantification

Verification

- Calibration tests using wafers from 5+ leading manufacturers
- Very good agreement with the same set of calibration data for all manufacturers
  - Independent of material properties (except thickness and offcut angle)
  - Recalibration / system check can be done based on XRT measurements, i.e., without additional KOH etching



Kranert et al., poster presentation @ECSCRM 2021





### **BPD** quantification – "The faster the better"

FastBPD approach using HyPix detector

- New version of XRTmicron available on request
  - Measurement times of ~5min per 150mm wafer
  - Identical results as for standard measurement
  - Insensitive to lattice curvature → wafers that fail on standard system can be measured
  - Very simple and fast alignment
  - For sufficiently tight specs → no alignment at all
  - Measurement reproducibility 5%
  - Tool2Tool: 5% or 50/cm<sup>2</sup>





### **BPD detection** SEMI Standard M93



SEMI M93 - Test Method for Quantifying Basal Plane Dislocation Density in 4H-SiC by X-Ray Diffraction Topography/Imaging

Member Price: \$138.00

Non-Member Price: \$180.00

Volume(s): Materials

Language: English

Type: Single Standards Download (.pdf)

SEMI Standards Copyright Policy/License Agreements

### **READY FOR INDUSTRY!**



### **XRT Toolbox - Defect analytics software**

- Experience from defect analysis condensed into XRT Toolbox
- Very fast defect analysis
   (<5 min for a full 150mm wafer scan)</li>
- Currently built for analysis of SiC wafers
- Can be adapted also for an analysis of defects in other materials
- Implemented evaluation routines fully compatible to SEMI Standards

### **READY FOR INDUSTRY!**





### **Summary**

- XRT has become a mature technique ready for application in industry
  - High accuracy, high reproducibility
  - SEMI-standardized procedures
  - High throughput
  - **Further developments under progress!**
  - Please be invited for further discussion!





### Center of Expertise for X-Ray-Topography



# ご清聴ありがとうございました。

## Thank you for your kind attention!

0





### Get in contact:

Dr.-Ing. Christian Reimann

**Rigaku Europe SE** 

Director Business Development christian.reimann@rigaku.com Dr. Assunta Vigliante **Rigaku Europe SE** Director Semiconductor Solutions assunta.vigliante@rigaku.com

0

**Dr. Christian Kranert** 

Fraunhofer IISB Group Manager Crystal- and Wafer Metrology christian.kranert@iisb.fraunhofer.de

1.1.1